



CARRIER CARDS FOR DXT-RAD RINGLETS

APPLICATION

The carrier card handles 4 pcs Thermo Scientific DXT-RAD ringlets and can be used in readers type Harshaw. The card is made of aluminum and the ringlets are fixed with 4 clamps. A semi automation for assembling and disassembling of the carrier card is available. The use of one of the Harshaw Standard TTPs (Publication DOSM-0-N-1201-001 dated 3. December 2002) is highly recommended.

SPECIFICATION

Material
Aluminum / Aluminum Alloy

Dimension
31 mm x 43,3mm x 1,5mm


Max. card temperature
90° Celsius / 20 min.

Reader acquisition temperature
max. 300° Celsius

Reader acquisition time
max. 20 sec.



Different colors available.



The card is provided with a serial / batch number.

CONTACT

Seibersdorf Labor GmbH
Radiation Protection Dosimetry
2444 Seibersdorf, Austria

Tel.: +43 50550 - 3000
Fax: +43 50550 - 3001

E-mail: dosimeterservice@seibersdorf-laboratories.at
Web: www.seibersdorf-laboratories.at/dosimeterservice